

▶ Mobility-Guard side rails

As part of our drive to promote patient autonomy, faster rehabilitation and reduce caregiver effort, the unique and innovative 3-position Mobility-Guard not only significantly facilitates patient self-repositioning and protects the patient from falls when fully raised, but also provides an intermediate position for safe patient egress.



Hill-Rom first designed split side rails in the 1950s in response to growing concerns regarding patient falls from beds. In 1974 we introduced bed controls embedded in the side rail. Hill-Rom introduced the first split side rails in Europe in 1995 with the Evolution bed. Since then, our pioneering TotalGuard® side rails have set a new standard in patient safety.

The side rail concept on the AvantGuard® 1600 bed is the latest in a long line of "firsts" developed by Hill-Rom.



▶ Making a difference to the lives of both patients and caregivers

The Hill-Rom AvantGuard® 1600 bed is designed to be Your Care Partner. Setting a new standard in performance to value, the AvantGuard® 1600 bed is ideally adapted to the demanding requirements of the care environment whilst enhancing the comfort and safety of patients and caregivers alike.

Because every facility is different, Hill-Rom proposes two models to meet your requirements:

AvantGuard® 1600 Ability Bed
Enhanced ergonomics and performance



AvantGuard® 1600 Mobility Bed
Enhanced patient mobilisation



designed by Twenty-Two, www.agntwentytwo.com

www.hill-rom.com

| | | | | | |
|----------------|----------------------|-------------------|---------------------|------------------|--|
| France | +33 (0)2 97 50 92 12 | Suisse/Schweiz | +44 (0)21/706 21 30 | Nordic Countries | +46 (0)8 564 353 60 |
| United Kingdom | +44 (0)1530 41 1000 | (deutschsprachig) | 706 21 38 | Export | contact your local Hill-Rom distributor |
| Deutschland | +49 (0)21 1/16450-0 | Österreich | +43 (0)2243/28550 | | or contact your Area manager via website |
| Nederland | +31 (0)347 32 35 32 | Ireland | +353 (0)1 413 6005 | | or call +1 812 934 8173 |
| Italia | +39 02-950541 | Iberia | +34 (0)93 6856000 | | |

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products. Manufactured by Hill-Rom SAS - Pluvigner, France © 2008 Hill-Rom Services, Inc. - ALL RIGHTS RESERVED.

USA & Canada
Western & Eastern Europe
Middle East / Africa / Indian Sub-Continent
Asia Pacific

Hill-Rom


Enhancing Outcomes for Patients and Their Caregivers.™

AvantGuard® 1600 Bed Your Care Partner



Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™



▶ Your Care Partner

Designed with a focus on ergonomics and safety, the Hill-Rom AvantGuard® 1600 bed delivers advanced capabilities that allow caregivers to meet the increasingly complex requirements inherent to acute care. State-of-the-art technologies and intelligent bed features contribute to more efficient and safer delivery of care, underscoring our commitment to "Clear Lungs, Safe Skin and No Falls™".

▶ Easy and Safe Patient Mobilisation

Mobilising the patient in the process of recovery is an essential part of care-giving. Caregiver injuries resulting from moving and handling of patients have been identified as one of the key ergonomics challenges within healthcare. With this in mind, the AvantGuard® 1600 bed features ergonomic capabilities designed to enhance the safety and efficiency of patient mobilisation.



▶ EasyChair® with 7 Bed Movements

At the touch of a button, the nurse can mobilise the patient to a fully seated position. Seven synchronised bed movements ensure optimal comfort and posture of the patient with minimal abdominal compression, shear and friction.

▶ Efficient Delivery of Care

The AvantGuard® 1600 bed offers "Nurse-Assist" features that ease the caregiver's workload and increase the quality of care: Intuitive and easy-to-use Point-of-Care® controls provide a single point of access to the bed's functions



The one-hand side rail release
Easily accessible for the caregiver to quickly lower the side rail.



A hip position indicator
allows for quick, easy and correct patient positioning.



Clocking movement
Ideally adapted to confined areas, the clocking movement of the side rails consumes no space at the bedside when raised or lowered.



User-friendly bed controls
allow the patient to reposition himself without assistance enhancing his autonomy and level of comfort.